

02/21/2024 LSRII/Fortessa Long Clean and PM Procedure – To be completed every 6 months.

Long Clean:

1. Add 0.5L bleach to sheath tank
2. Fill sheath tank to weld line with saline
3. Disconnect the 0.22uM fluid filter from sheath line
4. Directly connect sheath line (**Critical: Make sure filter is not connected!**) directly to sheath tank.
5. Connect sheath tank pressure line.
6. Open sheath line butterfly valve on side of instrument for 30secs to prime sheath line with 5% bleach saline.
7. Prime the system 5x
8. Place a tube of 10% bleach (bleach +DI) on to the system SIT.
9. Put system into a Low sample slow rate.
10. Press “Run” and allow the system to run for 60mins.
 - a. Keep an eye on the 10% bleach sample and refill as needed.

System Rinse

1. Dump remaining cleaning solution (5% bleach saline) from sheath tank.
2. Rinse sheath tank 5x with DI water from sink.
3. Refill sheath tank with saline.
4. Directly connect sheath line (san filter) directly to sheath tank.
5. Connect sheath tank pressure line.
6. Open sheath line butterfly valve on side of instrument for 30secs to rinse and prime sheath line with fresh saline.
7. Prime the system 5x
8. Place a tube of DI water onto the sit.
9. Put system into a Low sample slow rate.
10. Press “Run” and allow the system to run for 30mins.
 - a. Keep an eye on the Di water sample and refill as needed.

System Alignment

1. Adjust Blue, Red and Violet laser alignments to achieve <3.5% CV with CST Beads on primary detector.
2. Run CST Baseline (yearly)/Performance Check – check laser alignment, sample core alignment, detector health and fluidics.
3. **Save Performance Check report to log folder on desktop.**

Parts to replace

1. Sheath Filter
2. Waste Cap
3. Waste tubing
4. **Update Instrument log in log folder on desktop**